
SoloPath™ Balloon Expandable TransFemoral Introducer INSTRUCTIONS FOR USE

Description

Onset Medical Corporation's *SoloPath™* Balloon Expandable TransFemoral Introducer consists of a flexible, reinforced polymer sheath with a specially folded distal end (the Sheath) pre-mounted over a central balloon dilatation catheter (the Expander). The folded distal region of the Sheath is small in diameter, thus facilitating passage through the femoral artery. The *SoloPath* Assembly is inserted percutaneously into the femoral artery, over a guidewire, with the deflated Expander in place. The *SoloPath* is routed just past the iliac bifurcation into the aorta. The Expander balloon, when inflated with liquid, exerts controlled radial force, enlarging the folded distal region of the Sheath and surrounding anatomy. The Expander balloon is deflated and the Expander is removed leaving a large central lumen extending from the proximal end to the distal end of the Sheath. The Sheath is designed as a guide for catheters and/or devices introduced into the femoral iliac arteries.

Figure 1 – SoloPath™ With and Without Expander

Product Specifications

The Balloon Expandable TransFemoral Introducer is available with three (3) internal diameters and one (1) working length of 35cm. Refer to Table 1.

Table 1 – <i>SoloPath</i> [™] Balloon Expandable TransFemoral Introducer						
Model Number	Expanded Sheath I.D.	Expanded Sheath O.D.	Expandable Length	Dilator Tip Diameter	Working Length	Rated Inflation Pressure
	French	French	CM	French	CM	ATM/bar
STFI-1935	19	22	25	5.3	35	20
STFI-2035	20	23	25	5.3	35	20
STFI-2135	21	24	25	5.3	35	20

Warning: The *SoloPath*[™] Balloon Expandable TransFemoral Introducer is supplied sterile and is intended for single use only. Do not reuse, resterilize or reprocess this device. Do not use this device if product or sterile packaging is damaged.

STERILE EO

② Single Use Only

Device Construction

Intended for travel over a 0.038" or smaller guidewire.

A flexible Sheath composed of reinforced polymer is compressed, folded and wrapped around a central, removable balloon dilatation catheter (the Expander).

One radiopaque marker is located on the Sheath to indicate the end of the expandable portion of the Sheath.

The Expander balloon is non-compliant and size limited. The Expander has both an inflation lumen and a guidewire lumen.

The Expander is configured to be used with standard inflation devices having an appropriate pressure rating, a male Luer lock coupler, and a capacity of 25-CC.

The Sheath is terminated, at its proximal end, by a hub having a straight through hemostasis port capable of accepting catheters equal to the full sheath inner diameter.

An aspiration and infusion port, terminated with a 3-way stopcock, permits infusion of heparinized saline into the central port of the Sheath hub, prior to and during the procedure. The Expander incorporates a T connector, the through lumen port being terminated with a hemostasis valve for guidewire introduction, while the right angle port is for inflation of the Expander balloon. A separate aspiration and infusion side port, terminated with a 3-way stopcock, is also provided which permits fluid transfer into or from the guidewire lumen of the Expander.

Indications for Use

The *SoloPath*[™] Balloon Expandable TransFemoral Introducer is intended to be inserted percutaneously into the femoral artery, over a guidewire, and once expanded, to provide a guide for catheters and/or devices introduced into the femoral iliac arteries.

Caution: A thorough understanding of the technical principles, clinical applications and risks associated with cardiac transfemoral arterial access is required before attempting to use this device.

Contraindications

There are no known contraindications for the *SoloPath*[™] Balloon Expandable TransFemoral Introducer.

Precautions

- Inspect all components before use.
- Carefully read the instructions for use (IFU) before using this device to help reduce the potential dangers associated with the transfemoral technique such as air emboli and/or perforation of the iliac or aorta.
- Prior to inserting the device into the patient, insert the guidewire through the balloon Expander central lumen, checking and correcting for excessive resistance as the tip of the guidewire advances through the Sheath/Balloon Expander assembly.
- Frequently aspirate and flush, with heparinized saline, the Sheath, Expander, and any catheters to minimize the potential for air emboli and/or clot formation.
- Do not remove the balloon Expander or catheter from the Sheath rapidly or damage to the hemostasis valves may occur, causing blood leakage out through the valve.
- If resistance is met when advancing or withdrawing the guidewire or Sheath, determine the cause by fluoroscopy and correct before continuing with the procedure.

- Aspirate **only** from the Expander guidewire lumen or Sheath flush ports or Expander inflation port.
- Inject or provide saline flush only through the Expander or Sheath flush port or Expander inflation port.

Potential Complications

Complications associated with the use of the *SoloPath™* Balloon Expandable TransFemoral Introducer include but are not limited to, iliac perforation, aortic perforation, infection, post-procedural discomfort, air embolism, tissue trauma, hemorrhage, thromboembolism, thrombosis, injury to the vascular introduction site, and death.

Warnings

- Do not alter the device in any way.
- Federal Law (USA) restricts this device to sale by or on the order of a physician.
- Do not reuse this device. Thorough cleaning of biological and foreign material from the device is not possible. Adverse patient reactions may result from reuse of this device.
- Only those physicians who are specially trained in cardiac transfemoral arterial access procedures should use this device.
- Maintain continuous patient blood pressure monitoring throughout the procedure.
- Immediate availability of anterior-posterior and lateral fluoroscopy is required throughout the procedure.
- Always observe acceptable pressure tracings prior to advancing the dilator or any other component.
- Do not permit air to enter the Sheath. Remove components and make catheter exchanges slowly.
- Using the flushing ports, purge all air from the device prior to infusion.
- Provide a continuous drip of heparinized saline, under pressure, through the device while the Sheath remains within the vessel or cardiac structures.
- Fibrin may accumulate in or on the Sheath tip during the procedure. Aspirate following the removal of the Expander or catheters.

Suggested Procedure

Preparation

A standard preparation for endovascular access should be completed per hospital protocol. Standard fluoroscopic equipment should be available for use during the procedure. Proper radiological protection should be provided for all attending personnel.

Device Inspection

Using aseptic technique, remove the *SoloPath™* from its sterile pouch. Visually inspect the device to make sure there is no distortion or kinking in the shaft of the sheath or in the folded distal end and that a smooth taper exists between the distal end of the Sheath and the balloon Expander. The balloon Expander shaft is fixed into position within the folded Sheath body.

Caution: Do NOT attempt to separate the Sheath from the balloon Expander at any time prior to expansion of the Sheath within the iliac artery.

Device Preparation

Use aseptic technique for all steps of the procedure. Open the three-way stopcock on the Through Lumen Flush Port of the Sheath and flush with sterile, heparinized saline. Close the stopcock. Close the stopcocks to prevent saline loss or air embolism during the procedure.

Percutaneous Access

A standard Seldinger preparation and access are completed into the femoral artery using, for example, an 18-gauge thin wall access hollow needle. A 0.038" stiff guidewire with a floppy tip is advanced through the access needle and into the femoral artery. Route the 0.038" guidewire into the iliac artery, and just past the iliac bifurcation into the aorta under fluoroscopic guidance.

Device Insertion

Ensure that all ports have been flushed and primed with heparinized saline and are free of air. The *SoloPath™* with its integral balloon Expander is advanced as a single unit, under fluoroscopic control, over the 0.038" guidewire until its radiopaque marker is positioned just past the iliac bifurcation into the aorta. Inject radiopaque contrast media through the Expander flush port, as required, under fluoroscopic visualization to ensure correct placement.

Caution: The *SoloPath™* is to be advanced over the 0.038" guidewire as a single unit. Do not allow separation or misalignment of the central balloon Expander from the unexpanded Sheath prior to expansion within the vasculature.

Caution: Should resistance be encountered, cease advancing the *SoloPath™* until the cause of the resistance can be determined and corrected.

Caution: The balloon Expander is intended for inflation only after the unexpanded Sheath has been completely advanced to its final target location within the vasculature.

Caution: Do not attempt to refold an expanded sheath.

Sheath Expansion

A dilute 50% solution of radiopaque contrast media and sterile saline is prepared and approximately 15-CC are drawn up into a high-pressure balloon inflation syringe equipped with a pressure gauge (not supplied). Care should be taken to remove all air from the syringe, the Expander, and associated tubing. Attach the pressure line of the inflation device to the Expander balloon inflation port.

Under fluoroscopic control, the dilute contrast media is injected up to the maximum Rated Inflation Pressure per Table 1. It is normal to observe a drop in pressure as the Expander balloon progressively expands the folded distal section of the Sheath. The rated inflation pressure should be maintained for a minimum of 60 seconds to expand any "waist" that may remain along the length of the expanded Sheath.

Warning: Only use liquid to inflate the balloon Dilator. NEVER use gas, such as air, to inflate the balloon Dilator.

Deflated Expander Removal and Sheath Use

Apply suction, using the high-pressure syringe, to the inflation/deflation port on the balloon Dilator hub in order to deflate the balloon of the Expander. Remove the deflated Expander from the expanded Sheath slowly to minimize blood loss. Introduce the appropriate therapeutic device or catheter through the central working channel of the expanded Sheath.

Caution: Do NOT attempt to separate the balloon Expander from the Sheath while the Expander is under pressure.

Caution: Care should be taken to maintain the position of the expanded Sheath within the vasculature during the procedure. Advancing an already expanded Sheath must be avoided as this movement may not be possible and any attempt will potentially damage the vessel.

Sheath Removal

Remove any instrumentation from the Sheath slowly through the hemostasis valve, being careful to minimize blood loss. While maintaining hemostasis and using standard hospital procedure to control the percutaneous puncture site, gently withdraw and remove the expanded Sheath from the arterial circulation when the procedure has been completed.

How Supplied

Contents of unopened, undamaged package are sterile and non-pyrogenic.

STERILE EO

② Single Use Only

C€0297

This device is sterilized with ethylene oxide gas. It is intended for single use only. Do not use if the package is opened or damaged.

Storage

Store in a cool, dry place. Rotate inventory so that products are used prior to the sterilization expiration date on package label.

Warranty

Warranty and Limitations: Onset warrants to Customer that its products are free from defect in design, workmanship and materials. For its sterile products Onset further warrants to Customer that such products will remain sterile to the expiration date specified on the applicable product's label, provided that the original packaging remains intact, and is not compromised in any manner. This warranty shall not apply to any products, which have been re-sterilized, repaired, altered or modified in any manner, or to any products, which have been improperly stored, installed, operated, used or maintained in any manner. In the event of any breach of its warranty set forth above, Onset's sole obligation shall be to repair or replace, at its sole option, any product that Onset determines was defective in workmanship or materials at the time of shipment if notice thereof from Customer is received by Onset within the "use before" or expiration date, as applicable, described on such product's label. Except as expressly provided above, Customer assumes all responsibility and all other liability, whether based upon warranty, contract, negligence, or otherwise, for injury or damages resulting from handling, possession, use or misuse of any Onset product and agrees to defend and indemnify Onset for same. Because Onset has no control over the operation, use, or maintenance of its products, the selection of patients or their condition, and does not warrant the performance of its products, THE WARRANTY SET FORTH ABOVE IS EXPRESSLY ACKNOWLEDGED BY CUSTOMER TO BE IN LIEU OF ANY OTHER EXPRESS WARRANTY AND OF ANY OTHER OBLIGATION ON THE PART OF ONSET. The remedies set forth in this Warranty and Limitations shall be the exclusive remedy available to any person. No officer, employee, representative or any other agent of Onset has any authority to change or otherwise modify in any respect any of the foregoing, or to act in any manner to assume or bind Onset to any additional liability or responsibility in connection to its products. Customer's submission of any orders for any Onset product(s) shall be deemed irrevocable acceptance of all of the terms and conditions of the foregoing Warranty and Limitations.

EXCEPT FOR THE EXPRESS WARRANTY DESCRIBED ABOVE AND AS OTHERWISE SET FORTH HEREIN, ONSET GRANTS NO OTHER WARRANTIES OF ANY NATURE, DIRECTLY OR INDIRECTLY, EXPRESS OR IMPLIED, BY STATUTE OR OTHERWISE, REGARDING THE PRODUCTS, THEIR FITNESS OR SUITABILITY FOR ANY PURPOSE, THEIR QUALITY, CONDITION,

MERCHANTABILITY, OR OTHERWISE. ONSET'S LIABILITY UNDER THE WARRANTY SHALL BE LIMITED TO REPAIR OR REPLACEMENT OF PRODUCTS IN EXCHANGE FOR PRODUCTS RETURNED TO ONSET. IN NO EVENT SHALL ONSET BE LIABLE FOR ANY LOSS OF PROFITS, SPECIAL, CONSEQUENTIAL OR INCIDENTAL DAMAGES FOR BREACH OF LIABILITY.

U.S. and international patents pending
© Copyright 2009 Onset Medical Corporation
SoloPath™ is a trademark of Onset Medical Corporation

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Maximum Inflation Pressure

Manufacturer

Onset Medical Corporation
13900 Alton Parkway, Suite 120
Irvine, CA 92618 USA
949.716.1100
949.716.1101 Fax

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SoloPath™ ist eine Marke der Onset Medical Corporation