




Pathway™ Balloon Expandable Ureteral Access Sheath



The **Pathway™ Balloon Expandable Ureteral Access Sheath** is available in the following Dimension / Specifications:

Expanded Diameters	11Fr (ID) x 14Fr (OD) 12Fr (ID) x 15Fr (OD)	Model 255100, 255101, 255102 Model 255106, 255107, 255108
Guidewire Compatibility:	0.038" max	All Models
Working Length	28cm	Model 255100, 255106
	36cm	Model 255101, 255107
	46cm	Model 255102, 255108
Expandable Length	20cm	Model 255100, 255106
	25cm	Model 255101, 255102, 255107, 255108
Dilator Tip Dia.	4.4Fr	All Models
Expandable Dia. (Folded)	9.50Fr	Model 255100, 255101, 255102
	9.75Fr	Model 255106, 255107, 255108
Rated Dilator Pressure	25ATM	All Models
Hydrophilic Coating	Yes	All Models