



Pathway™ Balloon Expandable PCNL Sheath

The Next Generation for PCNL Procedures

A One-Step Percutaneous Nephrolithotomy Sheath

The procedure of choice for removal of large renal calculi is percutaneous nephrolithotomy (PCNL)

Advances in technique and equipment have allowed urologists to perform percutaneous stone removal with increasing efficacy and decreasing complications

A novel, single-step, balloon expandable, renal access device is now available, providing:

- Radial Dilation for Reduced Tissue Trauma
- General Access Including Upper Pole
- Accurate Guidewire Tracking and Guidewire Channel
- Easier Insertion
- Reduced Time
- Less Force to Insert
- Malleable Sheath Design May Allow for Easier Stone Removal

Onset Medical Corporation
 13900 Alton Parkway Suite 120
 Irvine, CA 92618 USA

Ordering Information

Email: cs@onsetmedical.com
 Phone: (949) 716-1100 Ext. 15
 Fax: (949) 716-1101

Pathway Access Sheath Catheter® is a Registered Trademark of Onset Medical Corporation

US and International Patents Pending

© Copyright 2008 Onset Medical Corporation



Ordering Information

Pathway™ Balloon Expandable PCNL Sheath							
Number	Description	Expanded Sheath ID		Expanded Sheath OD		Working Length	Overall Assembly Length
		French	mm	French	mm		
215304	Pathway PCNL 30Fr. x 17cm	30	10.0	33	11.0	17	23.4
215306	Pathway PCNL 30Fr. x 20cm	30	10.0	33	11.0	20	26.4
Number		Description					
C02219		Excel 25cc/30atm Sterile Inflation Device with Stopcock					

