




## **Pathway™ Balloon Expandable PCNL Sheath**



**The Pathway™ Balloon Expandable PCNL Sheath is available in the following Dimension / Specifications:**

Expanded Diameters	30Fr (ID) x 33Fr (OD)	Model 215304, 215306
Guidewire Compatibility:	0.038" max	All Models
Working Length	17cm	Model 215304
	20cm	Model 215306
Expandable Length	12cm	All Models
Expandable Dia.(Folded)	21Fr	All Models
Rated Dilator Pressure	25ATM	All Models
Hydrophilic coated tip	Yes	All Models