




## SoloPath™ Balloon Expandable TransFemoral Introducer



**The SoloPath™ Balloon Expandable TransFemoral Introducer is available in the following Dimension / Specifications:**

Expanded Diameters	19Fr (ID) x 22Fr (OD) 21Fr (ID) x 24Fr (OD)	Model 1935 Model 2135
Guidewire Compatibility:	0.38" max	All models
Working Length	35cm	All Models
Expandable Length	25cm	All Models
Dilator Tip Dia.	5.3Fr	All Models
Expandable Dia. (Folded)	≤ 15Fr	All Models
Rated Dilator Pressure	20ATM	All Models
Hydrophilic coated tip	Yes	All Models

