


SoloPath™ Balloon Expandable Transseptal Introducer



The SoloPath™ Balloon Expandable Transseptal Introducer is available in the following Dimension / Specifications:

Expanded Diameters	16Fr (ID) x 19Fr (OD) 18Fr (ID) x 21Fr (OD)	Model 1670 Model 1870
Guidewire Compatibility:	0.038" max	All Models
Working Length	70cm	All Models
Expandable Length	8.5cm	All Models
Dilator Tip Dia.	5.3Fr	All Models
Expandable Dia.(Folded)	≤13.5Fr	All Models
Rated Dilator Pressure	16ATM	All Models
Hydrophilic coated tip	Yes	All Models