



*Delivering the next generation of
"Intelligent" access sheath technology*


SoloPath™ Balloon Expandable Transseptal Introducer

**Controlled
Deployment
Technology:**

Provides the smallest sheath profile with the largest internal diameter of any sheath in its class.

Small insertion profile to protect against unnecessary trauma during insertion and septal crossing.

Enabling access for multiple and/or large device applications.



SoloPath™ Balloon Expandable Transseptal Introducer
Using Controlled Deployment Technology, the SoloPath™ Balloon Expandable Transseptal Introducer provides a low profile, single catheter access across the septum. Once expanded, the SoloPath provides a stable and secure conduit allowing large or multiple catheter access into the left atrium.

Onset Medical Corporation
13900 Alton Parkway, Suite 120
Irvine, CA 92618 USA

Ordering Information
Email: cs@onsetmedical.com
Phone: (949) 716-1100 Ext. 15
Fax: (949) 716-1101

SoloPath™ is a Registered Trademark of Onset Medical Corporation

US and International Patents Pending

© Copyright 2010 Onset Medical Corporation

SoloPath™ Balloon Expandable Transseptal Introducer						
Number	Description	I.D.	O.D.	Expandable Length	Tip Diameter	Working Length
SoloPath-1870	SoloPath Transseptal Introducer 18/21Fr x 70cm	18 Fr	21 Fr	8.5 cm	5.3 Fr	70 cm
SoloPath-1670	SoloPath Transseptal Introducer 16/19Fr x 70cm	16 Fr	19 Fr	8.5 cm	5.3 Fr	70 cm
Number		Description				
C02219		Excel 25cc/30atm Sterile Inflation Device with Stopcock				