




Expanding the Universe of Ureteral Access

Pathway™ Balloon Expandable Ureteral Access Sheath

**Controlled
Deployment
Technology:**

provides the smallest sheath profile with the largest internal diameter of any sheath in its class



- Radial dilation instead of axial shear.
- Reduced trauma during insertion.
- Reduced insertion force.
- Improved flow rates.
- Lower distal profile for crossing strictures.
- Increased flexibility with malleable sheath design.

Onset Medical Corporation
13900 Alton Parkway Suite 120
Irvine, CA 92618 USA

Ordering Information
Email: cs@onsetmedical.com
Phone: (949) 716-1100 Ext. 15
Fax: (949) 716-1101

Pathway Access Sheath Catheter® is a
Registered Trademark of Onset Medical
Corporation

US and International Patents Pending

© Copyright 2008 Onset Medical Corporation

Ordering Information

<i>Pathway™ Balloon Expandable Ureteral Access Sheath</i>						
Number	Description	I.D.	O.D.	Folded O.D.	Tip Diameter	Working Length
255100	Pathway Ureteral Sheath 11/14Fr x 28cm	11 Fr	14 Fr	9.5 Fr	4.4 Fr	28 cm
255101	Pathway Ureteral Sheath 11/14Fr x 36cm	11 Fr	14 Fr	9.5 Fr	4.4 Fr	36 cm
255102	Pathway Ureteral Sheath 11/14Fr x 46cm	11 Fr	14 Fr	9.5 Fr	4.4 Fr	46 cm
255106	Pathway Ureteral Sheath 12/15Fr x 28cm	12 Fr	15 Fr	9.75 Fr	4.4 Fr	28 cm
255107	Pathway Ureteral Sheath 12/15Fr x 36cm	12 Fr	15 Fr	9.75 Fr	4.4 Fr	36 cm
255108	Pathway Ureteral Sheath 12/15Fr x 46cm	12 Fr	15 Fr	9.75 Fr	4.4 Fr	46 cm
Number		Description				
C02219		Excel 25cc/30atm Sterile Inflation Device with Stopcock				